

**Ethox** TUM-E-VAC® GASTRIC LAVAGE KIT

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK



- Single use lavage kit
- Non-sterile and latex free

**Lavage Tubing**  
Radiopaque strip confirms proper tube placement



**48" Tubing Length**  
Reference marks at 18", 20" and 22" from distal end



**Optional charcoal**  
See number 11



**INSTRUCTIONS FOR USE**

1. Close both clamps on lavage kit.
2. Fill or spike irrigation solution bag.
3. Hang solution bag above patient.
4. Position drainage bag below patient, confirm clamp and outlet bag are closed.
5. Insert lavage tube into patient using approved standard hospital procedure.
  - 5.1 Using standard X-ray technique, confirm proper placement of lavage tube.
6. Fill the irrigation tube by opening clamp located between solution bag and connector. Close clamp when tube is full.
7. Cut universal connector to desired length and diameter (Please reference tag on drainage tube for cutting universal connector to proper length). Attach to lavage tube. (Connector will fit lavage tube sizes 18 Fr. to 40 Fr.) Keep connector below patient's stomach level.
8. To start lavage procedure, open irrigation tube clamp until desired volume of fluid has been administered. Close clamp.
9. To drain stomach fluid, open drain tube clamp. When drainage is complete, close clamp.
  - 9.1 If fluid is not draining properly, the lavage tube may be obstructed.
  - 9.2 To clear lavage tube, disconnect tube at connector. Using the provided piston syringe, inject irrigation solution into lavage tube. Attach lavage tube to connector.
  - 9.3 If lavage tube remains obstructed, remove tube from patient, correct situation, and repeat step 5.
10. Repeat steps 8 and 9 as required.

**Model Numbers 2057, 2058, 2059, 2081 only:**

11. After lavage is complete, snip tip on charcoal tube, remove red cap on connector and thread charcoal tube into connector. Slide Twist-Lok key onto the end of the charcoal tube and twist to dispense charcoal. After charcoal is dispensed, slide the Twist-Lok to the side to lock under the tube.

NOTE: The enclosed syringe is an accessory for charcoal delivery, lavage tube placement, or clearing lavage tube blockage. It has no measuring function.