

Ethox

INFU-SURG® PRESSURE INFUSION BAG

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

The latex free Infu-Surg® Pressure Infusion Bag is a safe, high-quality, reliable device for intravenous and intra-arterial infusion treatment, including A-line pressure monitoring. The Infu-Surg is a cost-effective, trusted product that medical professionals have relied on for decades. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use

INFU-SURG® PRESSURE INFUSION BAG

- Pressure relief valve prevents over inflation (330mmHg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300mmHg)
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested

CAT#	VOLUME	DESCRIPTION	PK
4005H	500ml	Infu-Surg® with Hook	5
4010H	1000ml	Infu-Surg® with Hook	5
4030	3000ml	Infu-Surg®	5

