

SURGICAL

Ethox

SURGI-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Surgi-Kit® Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. This disposable tool maintains the sterile operating field and its simplicity eliminates cleaning procedures. Available sterile and non-sterile, choose from a variety of sizes and pocket configurations.



Keep instruments within arm's reach by attaching instrument pocket to operating drape; select sizes available with peel-n-stick tape.

SURGI-KIT® INSTRUMENT AND SUTURE HOLDER


- Keeps instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Available with tape to adhere to operating drape
- Available sterile and non-sterile

STERILE


CAT#	POCKET SIZE	DESCRIPTION	PK
SK-100	19cm x 23cm	8 pockets	50
SK-150	19cm x 11cm	4 pockets	100
SK-200	28cm x 18cm	2 pockets	50
SK-300	19cm x 23cm	8 pockets + tape	50
SK-350	19cm x 11cm	4 pockets + tape	75
SK-400	28cm x 18cm	2 pockets + tape	50
SK-650	5cm x 19cm	10 pockets + tape	50

NON-STERILE


CAT#	POCKET SIZE	DESCRIPTION	PK
SK-100N	19cm x 23cm	8 pockets	500
SK-150N	19cm x 11cm	4 pockets	500
SK-200N	28cm x 18cm	2 pockets	500
SK-300N	19cm x 23cm	8 pockets + tape	500
SK-350N	19cm x 11cm	4 pockets + tape	500
SK-400N	28cm x 18cm	2 pockets + tape	500
SK-650N	5cm x 19cm	10 pockets + tape	500




8 Pocket
SK-100
SK-100N




4 Pocket
SK-150
SK-150N




2 Pocket
SK-200
SK-200N



8 Pocket + Tape
SK-300
SK-300N



4 Pocket + Tape
SK-350
SK-350N



2 Pocket + Tape
SK-400
SK-400N



10 Pocket + Tape
SK-650
SK-650N





LAP-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Lap-Kit® Laparoscopic Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. Available sterile and non-sterile, choose from a variety of sizes and pocket configurations.



LAP-KIT® LAPAROSCOPIC INSTRUMENT HOLDER


- Keeps laparoscopic instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Available with tape to adhere to operating drape
- Available sterile and non-sterile

STERILE


CAT#	POCKET SIZE	DESCRIPTION	PK
LK-700	6.5cm x 30.5cm	2 Pockets + Tape	25
LK-800	3.2cm x 30.5cm	4 Pockets + Tape	25
LK-900	3.2cm x 30.5cm (2), 6.5cm x 30.5cm	3 Pockets + Tape	25

NON-STERILE


CAT#	POCKET SIZE	DESCRIPTION	PK
LK-700N	6.5cm x 30.5cm	2 Pockets + Tape	200
LK-900N	3.2cm x 30.5cm (2), 6.5cm x 30.5cm	3 Pockets + Tape	200



2 Pockets
LK-700
LK-700N



4 Pockets
LK-800



3 Pockets
LK-900
LK-900N