

SunMed FIBER OPTIC GREENLINE®/D™
ALL-METAL BLADE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

Fiber Optic All-Metal laryngoscope blades solve contamination problems and eliminate the cost and time spent cleaning blades and returning them for use. These surgical grade 303/304 stainless steel disposables also answer the professional's request for a non-plastic disposable blade. The flexible fiber optic bundle is protected in a black plastic sheath. All sizes are compatible with any reusable "Green System" Fiber Optic handle and comply with ISO 7376 Green Standard.



ALL-METAL MILLER LARYNGOSCOPE BLADES

CAT#	SIZE	DESCRIPTION	A	B	PK
5-5336-20	00	Premature	67mm	10mm	20
5-5336-00	0	Neonate	77mm	10mm	20
5-5336-01	1	Infant	92mm	10mm	20
5-5336-02	2	Child	100mm	12mm	20
5-5336-03	3	Medium Adult	130mm	14mm	20
5-5336-04	4	Large Adult	155mm	18mm	20

ALL-METAL MACINTOSH LARYNGOSCOPE BLADES

CAT#	SIZE	DESCRIPTION	B	A	PK
5-5335-01	1	Infant	92mm	22mm	20
5-5335-02	2	Child	100mm	21mm	20
5-5335-03	3	Medium Adult	130mm	22mm	20
5-5335-04	4	Large Adult	155mm	25mm	20

