

AIRWAY MANAGEMENT



ENDOTRACHEAL TUBE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE | CE MARK



15mm adapter with wings and knurling facilitates connector removal when rocking adapter loose

Tapered ID inside connector reduces flow turbulence into tube



Textured Aluminum, low friction stylet slides freely in and out of tube

Ultra-secure adapter helps prevent disconnection

Gas tight female luer slip inflation valve

Radiopaque embedded strip enables tube position detection with x-ray

Sensitive pilot balloon indicates status of cuff

Bold CM depth markings are easy to read in low light conditions

Polished 45° beveled tip helps reduce trauma on insertion

Polished Murphy Eye facilitates ventilation if tip becomes occluded

No ridge where the cuff meets the tube allows for more gentle insertion and removal at vocal cords

Crystal clear PVC allows immediate detection of foreign material in tube

Clearly visible mm OD and ID size marking helps proper tube size selection

Two positioning marks aid in correct tube depth

Barrel cuff with generous surface area lowers pressure in trachea capillaries while minimizing folds in cuff

Low pressure, high volume cuff provides a secure, comfortable seal in trachea



SunMed ENDOTRACHEAL TUBE

SunMed's endotracheal tubes are designed for either oral or nasal use with a smooth Murphy Eye as an alternative ventilation port. Made of latex free, flexible PVC, that are available cuffed, uncuffed, with or without stylet. Cuffed tubes feature a high volume, low pressure barrel cuff. A radiopaque strip is embedded to assess exact location of tube. All have a standard 15mm ISO connector. These airways are individually packaged.



ENDOTRACHEAL TUBE - WITH STYLET

CAT#	ID	SIZE	LENGTH	BOX
1-7343-30	3.0mm	12Fr	165mm	10
1-7343-35	3.5mm	14Fr	185mm	10
1-7343-40	4.0mm	16Fr	210mm	10
1-7343-45	4.5mm	18Fr	225mm	10
1-7343-50	5.0mm	20Fr	245mm	10
1-7343-55	5.5mm	22Fr	275mm	10
1-7343-60	6.0mm	24Fr	285mm	10
1-7343-65	6.5mm	26Fr	295mm	10
1-7343-70	7.0mm	28Fr	310mm	10
1-7343-75	7.5mm	30Fr	315mm	10
1-7343-80	8.0mm	32Fr	330mm	10
1-7343-85	8.5mm	34Fr	330mm	10
1-7343-90	9.0mm	36Fr	330mm	10
1-7343-95	9.5mm	38Fr	330mm	10
1-7343-10	10.0mm	40Fr	330mm	10

ENDOTRACHEAL TUBE - CUFFED

CAT#	ID	SIZE	LENGTH	BOX
1-7333-30	3.0mm	12Fr	165mm	10
1-7333-35	3.5mm	14Fr	185mm	10
1-7333-40	4.0mm	16Fr	210mm	10
1-7333-45	4.5mm	18Fr	225mm	10
1-7333-50	5.0mm	20Fr	245mm	10
1-7333-55	5.5mm	22Fr	275mm	10
1-7333-60	6.0mm	24Fr	285mm	10
1-7333-65	6.5mm	26Fr	295mm	10
1-7333-70	7.0mm	28Fr	310mm	10
1-7333-75	7.5mm	30Fr	315mm	10
1-7333-80	8.0mm	32Fr	330mm	10
1-7333-85	8.5mm	34Fr	330mm	10
1-7333-90	9.0mm	36Fr	330mm	10
1-7333-95	9.5mm	38Fr	330mm	10
1-7333-10	10.0mm	40Fr	330mm	10



CUFF INFLATION BULB

- Quickly inflates and deflates endotracheal tube cuffs
- Safe - no syringes
- One-way check valve
- Pocket size - one-hand operation
- Latex free, non-sterile

CAT#	DESCRIPTION	BOX
4-0101-00	Cuff Inflation Bulb	10



ENDOTRACHEAL TUBE - UNCUFFED

CAT#	ID	SIZE	LENGTH	BOX
1-7330-20	2.0mm	8Fr	140mm	10
1-7330-25	2.5mm	10Fr	145mm	10
1-7330-30	3.0mm	12Fr	165mm	10
1-7330-35	3.5mm	14Fr	185mm	10
1-7330-40	4.0mm	16Fr	210mm	10
1-7330-45	4.5mm	18Fr	225mm	10
1-7330-50	5.0mm	20Fr	245mm	10
1-7330-55	5.5mm	22Fr	275mm	10

