


SunMed FIBER OPTIC GREENLINE®
GERMAN PROFILE BLADE

LATEX FREE | CE MARK

Fiber Optic German Profile laryngoscope blades are built with an integral, extra large fiber optic bundle providing a more cool, white light. This integral design enables superior cleaning by eliminating cracks and crevices while preventing cross contamination. These blades are made of 303/304 surgical stainless steel, not chrome plated.



MACINTOSH GERMAN PROFILE

CAT#	SIZE	DESCRIPTION	A	B	PK
5-5272-01	1	Infant	93mm	15mm	1
5-5272-02	2	Child	114mm	20mm	1
5-5272-13	3	Strong Curve Medium Adult	135mm	20mm	
5-5272-03	3	Medium Adult	135mm	20mm	1
5-5272-04	4	Large Adult	155mm	20mm	1

MILLER GERMAN PROFILE

CAT#	SIZE	DESCRIPTION	A	B	PK
5-5273-20	00	Premature	65mm	11mm	1
5-5273-00	0	Newborn	77mm	12mm	1
5-5273-01	1	Infant	110mm	12mm	1
5-5273-15	1.5	Small Child	130mm	14mm	1
5-5273-02	2	Child	152mm	14mm	1
5-5273-03	3	Medium Adult	192mm	15mm	1
5-5273-04	4	Large Adult	207mm	16mm	1

