


ANESTHESIA

SunMed FIBER OPTIC GREENLINE® / D™ BLADE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

SunMed's Fiber Optic GreenLine®/D™ disposable laryngoscope blades help solve contamination problems and eliminate the cost and time spent cleaning blades and returning them for use. These stainless steel blades also answer the professional's request for a non-plastic disposable blade. The polished acrylic stem produces exceptional illumination. The extended safety heel prevents blade-to-handle contact. All sizes and styles are compatible with Green System fiber optic handles and are ISO 7376 compliant.



MILLER GREENLINE® / D™

CAT#	SIZE	DESCRIPTION	A	B	PK
5-5333-20	00	Premature	67mm	10mm	20
5-5333-00	0	Neonate	77mm	10mm	20
5-5333-01	1	Infant	102mm	10mm	20
5-5333-02	2	Child	153mm	12mm	20
5-5333-03	3	Medium Adult	195mm	14mm	20
5-5333-04	4	Large Adult	205mm	18mm	20

MACINTOSH GREENLINE® / D™

CAT#	SIZE	DESCRIPTION	A	B	PK
5-5332-01	1	Infant	92mm	22mm	20
5-5332-02	2	Child	100mm	21mm	20
5-5332-03	3	Medium Adult	130mm	22mm	20
5-5332-35	3.5	Ext. Medium Adult	144mm	22mm	20
5-5332-04	4	Large Adult	155mm	25mm	20

