

SunMed FIBER OPTIC GREENLINE® / D™ HANDLE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

GreenLine®/D™ disposable LED laryngoscope handles are ideal for use with GreenLine® Fiber Optic laryngoscope blades. This single use handle helps address secondary infection concerns and costs associated with reprocessing.



5-0236-69

5-0236-70

5-5338-11

5-5338-10

5-5338-09

5-5338-12

GREENLINE® / D™ DISPOSABLE LED HANDLE

- One complete, fully-contained unit
- Ideal for use with GreenLine®/D™ disposable blades
- LED illumination provides exceptional, clear viewing
- Complies with ISO 7376 Standard
- Batteries included

CAT#	SIZE	DESCRIPTION	STERILE	PK
5-0236-69	Large	Ribbed Grip	Clean Pkg	20
5-0236-70	Penlite	Smooth Grip	Clean Pkg	20
5-5338-09	Medium	Textured Grip	Yes	20
5-5338-10	Penlite	Textured Grip	Yes	20
5-5338-11	Stubby	Textured Grip	Yes	20
5-5338-12	Standard	Textured Grip	Yes	20

