

SunMed IGuard™ EYE PROTECTOR

LATEX FREE | STERILE | SINGLE USE | CE MARK


The IGuard provides comfortable, effective patient eye protection while in the ER, surgery and during general procedures. Available in single or double thickness foam, the IGuard stays in place with non-allergenic adhesive.




IGUARD EYE PROTECTOR

- Rigid, clear protective patient eye cover
- Non-allergenic self-adhesive foam cushion for fast and accurate application
- Thick foam for patient comfort even in prone position
- Improved one piece backing with "courtesy tabs"

CAT#	DESCRIPTION	PK
9-0210-00	Adult Eye Protector	25
9-0210-21	Pediatric Eye Protector	25
9-0210-20	Double Foam Eye Protector	25



9-0210-20



9-0210-00