

AIRWAY MANAGEMENT

SunMed ORAL AIRWAY

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Latex free oral Guedel airways are designed with built-in bite blocks. Both styles are color coded for accurate and quick size identification. These airways are non-sterile and individually packaged; sets are ideal for crash carts.



COLOR CODED GUEDEL AIRWAY

- Firm, polyethylene plastic

CAT#	SIZE	DESCRIPTION	COLOR	PK
1-1504-40	4cm	Newborn, 40mm	Pink	10
1-1504-50	5cm	Infant, 50mm	Blue	10
1-1504-60	6cm	Small Child, 60mm	Black	10
1-1504-70	7cm	Child, 70mm	White	10
1-1504-80	8cm	Small Adult, 80mm	Green	10
1-1504-90	9cm	Med. Adult, 90mm	Yellow	10
1-1504-99	10cm	Large Adult, 100mm	Red	10
1-1504-11	11cm	X-Large Adult, 110mm	Orange	10
1-1505-40	4cm	Newborn, 40mm	Pink	50
1-1505-50	5cm	Infant, 50mm	Blue	50
1-1505-60	6cm	Small Child, 60mm	Black	50
1-1505-70	7cm	Child, 70mm	White	50
1-1505-80	8cm	Small Adult, 80mm	Green	50
1-1505-90	9cm	Med. Adult, 90mm	Yellow	50
1-1505-99	10cm	Large Adult, 100mm	Red	50
1-1505-11	11cm	X-Large Adult, 110mm	Orange	50

SETS

CAT#	SIZES INCLUDED	PK
1-1504-00	40, 50, 60, 70, 80, 90, 100, 110mm	8/pack

COLOR CODED SOFT PVC GUEDEL AIRWAY

- Soft, clear PVC plastic

CAT#	SIZE	DESCRIPTION	COLOR	PK
1-5010-30	3cm	Premature, 30mm	Clear	10
1-5010-20	4cm	Newborn, 40mm	Pink	10
1-5010-10	5cm	Infant, 50mm	Blue	10
1-5010-01	6cm	Small Child, 60mm	Black	10
1-5010-02	7cm	Child, 70mm	White	10
1-5010-03	8cm	Small Adult, 80mm	Green	10
1-5010-04	9cm	Med. Adult, 90mm	Yellow	10
1-5010-05	10cm	Large Adult, 100mm	Red	10
1-5010-06	11cm	X-Large Adult, 110mm	Orange	10
1-1521-30	3cm	Premature, 30mm	Clear	50
1-1521-20	4cm	Newborn, 40mm	Pink	50
1-1521-10	5cm	Infant, 50mm	Blue	50
1-1521-01	6cm	Small Child, 60mm	Black	50
1-1521-02	7cm	Child, 70mm	White	50
1-1521-03	8cm	Small Adult, 80mm	Green	50
1-1521-04	9cm	Med. Adult, 90mm	Yellow	50
1-1521-05	10cm	Large Adult, 100mm	Red	50
1-1521-06	11cm	X-Large Adult, 110mm	Orange	50

SETS

CAT#	SIZES INCLUDED	PK
1-5010-98	30, 40, 50, 60, 70, 80, 90, 100, 110mm	9/pack



SunMed ORAL AIRWAY

SunMed's latex free oral airways are made of firm plastic with built-in bite blocks. The Berman airways are color coded for accurate and quick size identification. These single use airways are non-sterile; sets are ideal for crash carts.



COLOR CODED BERMAN AIRWAY

- Firm, polyethylene plastic
- Vented septum

CAT#	SIZE	DESCRIPTION	COLOR	PK
1-1508-40	4cm	Newborn, 40mm	Pink	10
1-1508-50	5cm	Infant, 50mm	Blue	10
1-1508-60	6cm	Small Child, 60mm	Black	10
1-1508-70	7cm	Child, 70mm	White	10
1-1508-80	8cm	Small Adult, 80mm	Green	10
1-1508-90	9cm	Med. Adult, 90mm	Yellow	10
1-1508-99	10cm	Large Adult, 100mm	Red	10
1-1508-11	11cm	X-Large Adult, 110mm	Orange	10

1-1506-40	4cm	Newborn, 40mm	Pink	50
1-1506-50	5cm	Infant, 50mm	Blue	50
1-1506-60	6cm	Small Child, 60mm	Black	50
1-1506-70	7cm	Child, 70mm	White	50
1-1506-80	8cm	Small Adult, 80mm	Green	50
1-1506-90	9cm	Med. Adult, 90mm	Yellow	50
1-1506-99	10cm	Large Adult, 100mm	Red	50
1-1506-11	11cm	X-Large Adult, 110mm	Orange	50

SETS

CAT#	SIZES INCLUDED	PK
1-1508-00	40, 50, 60, 70, 80, 90, 100, 110mm	10 sets



TRADITIONAL GUEDEL AIRWAY

- Firm, polyethylene plastic
- One piece construction
- Built-in bite block
- Bulk packaged

CAT#	SIZE	DESCRIPTION	COLOR	PK
1-1501-50	5cm	Infant, 50mm	Blue	10
1-1501-60	6cm	Small Child, 60mm	Black	10
1-1501-70	7cm	Child, 70mm	White	10
1-1501-80	8cm	Small Adult, 80mm	Green	10
1-1501-90	9cm	Med. Adult, 90mm	Yellow	10
1-1501-99	10cm	Large Adult, 100mm	Red	10

1-1500-50	5cm	Infant, 50mm	Blue	50
1-1500-60	6cm	Small Child, 60mm	Black	50
1-1500-70	7cm	Child, 70mm	White	50
1-1500-80	8cm	Small Adult, 80mm	Green	50
1-1500-90	9cm	Med. Adult, 90mm	Yellow	50
1-1500-99	10cm	Large Adult, 100mm	Red	50



WILLIAMS AIRWAY INTUBATOR

- Guides fiber optic bronchoscopy placement
- Built-in bite block
- Bulk package

CAT#	SIZE	DESCRIPTION	PK
1-1503-09	9cm	Adult Female	10
1-1503-10	10cm	Adult Male	10

