

EMERGENCY


SunMed

CATHETER

STERILE | SINGLE USE | DISPOSABLE | CE MARK


Foley pediatric and adult catheters are designed with a smooth tapered tip for easy insertion into the urethra. They also feature symmetrical balloons that expand equally in all directions and valves that permit simple inflation and deflation.

Available in latex silicone coated or 100% silicone.



7-6503 Series

7-6505 Series



7-6515 Series

7-6513 Series

LATEX SILICONE COATED FOLEY CATHETER

- Smooth tapered tip facilitates easy insertion into urethra
- Drainage eyes are accurately formed to permit effective drainage
- Symmetrical balloon expands equally in all directions and efficiently retains bladder
- Color-coded sleeves for easy and rapid size identification
- Plastic valve permits simple inflation and deflation

CAT#	SIZE	DESCRIPTION	PK
7-6503-08	8Fr	Pediatric, 2 ways, 3cc	10
7-6503-10	10Fr	Pediatric, 2 ways, 3cc	10
7-6505-12	12Fr	Adult, 2 ways, 5cc	10
7-6505-14	14Fr	Adult, 2 ways, 5cc	10
7-6505-16	16Fr	Adult, 2 ways, 5cc	10
7-6505-18	18Fr	Adult, 2 ways, 5cc	10
7-6505-20	20Fr	Adult, 2 ways, 5cc	10
7-6505-22	22Fr	Adult, 2 ways, 5cc	10

100% SILICONE FOLEY CATHETER (LATEX FREE)

- Silicone – Inert material preferred for medical devices, less urethral irritation and stricture versus latex material
- Less encrustation compared to latex catheters; 100% latex-free
- Smooth tapered tip facilitates easy insertion into urethra
- Symmetrical balloon expands equally in all directions and efficiently retains bladder
- Color-coded sleeves for easy and rapid size identification
- Plastic Bspak spring-loaded valve permits simple inflation and deflation
- More efficient than rubber valves

CAT#	SIZE	DESCRIPTION	PK
7-6513-08	8Fr	Pediatric, 2 ways, 3cc	10
7-6513-10	10Fr	Pediatric, 2 ways, 3cc	10
7-6515-12	12Fr	Adult, 2 ways, 5cc	10
7-6515-14	14Fr	Adult, 2 ways, 5cc	10
7-6515-16	16Fr	Adult, 2 ways, 5cc	10
7-6515-18	18Fr	Adult, 2 ways, 5cc	10
7-6515-20	20Fr	Adult, 2 ways, 5cc	10
7-6515-22	22Fr	Adult, 2 ways, 5cc	10

