




CONVENTIONAL MACINTOSH LARYNGOSCOPE BLADE

LATEX FREE | CE MARK

SunMed's English profile conventional laryngoscopes are available in both Macintosh and Miller designs. Constructed of 303/304 surgical stainless steel, these blades come with either LED or clear lamps. The LED lamp emits a cool white illumination. This durable performance plus low current use increases battery life.



MACINTOSH BLADES WITH LED LAMP

- Macintosh laryngoscope design is predominant choice among curved blades
- Flange extends all the way down to distal tip
- Soft matte finish virtually eliminates reflection and glare
- Cool white LED illumination delivers 35,000 hours of use; durable performance plus low current use increases battery life

ENGLISH PROFILE


CAT#	SIZE	DESCRIPTION	A	B
5-5051-00	0	Neonate	80mm	17mm
5-5051-01	1	Infant	92mm	14mm
5-5051-02	2	Child	100mm	17mm
5-5051-03	3	Medium Adult	130mm	18mm
5-5051-04	4	Large Adult	155mm	18mm
5-5051-05	5	X-Large Adult	175mm	24mm

AMERICAN PROFILE

CAT#	SIZE	DESCRIPTION	A	B
5-5052-00	0	Neonate	80mm	17mm
5-5052-01	1	Infant	92mm	14mm
5-5052-02	2	Child	100mm	17mm
5-5052-03	3	Medium Adult	130mm	18mm
5-5052-35	3.5	Ext. Medium Adult	144mm	22mm
5-5052-04	4	Large Adult	159mm	22mm

LEFT-HANDED BLADE

CAT#	SIZE	DESCRIPTION	A	B
5-5053-03	3	Medium Adult Left	130mm	18mm



SunMed CONVENTIONAL MILLER LARYNGOSCOPE BLADE

LATEX FREE | CE MARK

SunMed’s English profile conventional laryngoscopes are available in both Macintosh and Miller designs. Constructed of 303/304 surgical stainless steel, these blades come with either LED or clear lamps. The LED lamp emits a cool white illumination. This durable performance plus low current use increases battery life.



MILLER BLADES WITH LED LAMP

- Miller laryngoscope design is predominant choice among straight blades
- Satin finish virtually eliminates reflection and glare
- Cool white LED illumination delivers 35,000 hours of use; durable performance plus low current use increases battery life

CLEAR HALOGEN LAMP

CAT#	SIZE	DESCRIPTION	A	B
5-5061-00	0	Neonate	77mm	9mm
5-5061-01	1	Infant	102mm	9mm
5-5061-02	2	Child	153mm	11mm
5-5061-03	3	Medium Adult	195mm	11mm
5-5061-04	4	Large Adult	205mm	16mm

LED LAMP

CAT#	SIZE	DESCRIPTION	A	B
5-5861-00	0	Neonate	77mm	9mm
5-5861-01	1	Infant	102mm	9mm
5-5861-02	2	Child	153mm	11mm
5-5861-03	3	Medium Adult	195mm	11mm
5-5861-04	4	Large Adult	205mm	16mm

