

ANESTHESIA

**SunMed** FIBER OPTIC GREENLINE® HANDLE

LATEX FREE | CE MARK

Choose from a variety of handles for use with SunMed’s fiber optic laryngoscope reusable blades. All are compatible with Green System blades and comply with ISO 7376 Green Standard.



**GREENLINE® HANDLE**

- Chrome plated brass and surgical stainless steel handles are compatible with all “Green System” fiber optic blades.
- Bright Halogenated Xenon lamp

**CHROME PLATED**

CAT#	SIZE	BATTERIES*	BOX
5-0236-09	Medium Handle	2 “C”	1
5-0236-10	Penlite Handle	2 “AA”	1
5-0236-11	Stubby Handle	2 “AA”	1
5-0236-15	Large Handle	2 “D”	1

**SURGICAL STAINLESS STEEL**

CAT#	SIZE	BATTERIES*	BOX
5-0236-39	Medium Handle	2 “C”	1
5-0236-30	Penlite Handle	2 “AA”	1
5-0236-31	Stubby Handle	2 “AA”	1
5-0236-35	Large Handle	2 “D”	1
5-0236-38	Mini Handle	Lithium “123”	1
5-0236-36	Micro Handle	Lithium “123”	1

**PREMIUM LED HANDLE**

- Premium LED handles have a lifetime warranty
- Low power use prolongs battery life
- Compliant with ISO standard 7376



CAT#	SIZE	BATTERIES	BOX
5-0236-80	Premium LED Penlite	2 “AA”	1
5-0236-89	Premium LED Medium	2 “AA”	1

**SUNBRITE™ LED HANDLE**

- Eliminates any heat in handle
- Cool, white illumination
- Low power use prolongs battery life



CAT#	SIZE	BATTERIES*	BOX
5-0236-90	Penlite Handle	2 “AA”	1
5-0236-91	Stubby Handle	2 “AA”	1
5-0236-99	Medium Handle	2 “C”	1

\*Batteries not included

