



CONTRACT NUMBER MS4810

September 1, 2017 - September 30, 2019



Available with
economical **mesh** panel
or **clear** panel



Ethox

PRESSURE INFUSION BAGS


- Pressure relief valve prevents over inflation (330mmHg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300mmHg)
- Incredibly reliable – 100% tested



Clear solution bag panel enables patient label and barcode scanning without having to remove bag.



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand.



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation.

Novaplus™ Pressure Infusion Bag

Mesh Panel		Clear Panel	
Series	Size	Series	Size
V4005H	500cc	V4050	500cc
V4010H	1000cc	V4100	1000cc
4030*	3000cc		

*Ethox brand pressure infusion bag

Product details available
at www.Sun-Med.com