

# SURGI-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK


Surgi-Kit® Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. This disposable tool maintains the sterile operating field and its simplicity eliminates cleaning procedures. Available sterile and non-sterile, choose from a variety of sizes and pocket configurations.




Keep instruments within arm's reach by attaching instrument pocket to operating drape; select sizes available with peel-n-stick tape.

## SURGI-KIT® INSTRUMENT AND SUTURE HOLDER


- Keeps instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Available with tape to adhere to operating drape
- Available sterile and non-sterile




8 Pocket  
SK-100  
SK-100N




4 Pocket  
SK-150  
SK-150N




2 Pocket  
SK-200  
SK-200N




8 Pocket + Tape  
SK-300  
SK-300N



4 Pocket + Tape  
SK-350  
SK-350N



2 Pocket + Tape  
SK-400  
SK-400N



10 Pocket + Tape  
SK-650  
SK-650N

## STERILE

ITEM	DESCRIPTION	SIZE	PK
SK-100	8 pockets	19 cm x 23 cm	50
SK-150	4 pockets	19 cm x 11 cm	100
SK-200	2 pockets	28 cm x 18 cm	50
SK-300	8 pockets + tape	19 cm x 23 cm	50
SK-350	4 pockets + tape	19 cm x 11 cm	75
SK-400	2 pockets + tape	28 cm x 18 cm	50
SK-650	10 pockets + tape	5 cm x 19 cm	50

## NON-STERILE

ITEM	DESCRIPTION	SIZE	PK
SK-100N	8 Pockets	19 cm x 23 cm	50
SK-150N	4 Pockets	19 cm x 11 cm	100
SK-200N	2 Pockets	28 cm x 18 cm	50
SK-300N	8 Pockets + Tape	19 cm x 23 cm	50
SK-350N	4 Pockets + Tape	19 cm x 11 cm	75
SK-400N	2 Pockets + Tape	28 cm x 18 cm	50
SK-650N	10 Pockets + Tape	5 cm x 19 cm	50

# LAP-KIT®

LATEX FREE | SINGLE USE | DISPOSABLE | CE MARK

Lap-Kit® Laparoscopic Instrument Holders are disposable pockets for sterile delivery, common in surgical applications. Available sterile and non-sterile, choose from a variety of sizes and pocket configurations.



## LAP-KIT® LAPAROSCOPIC INSTRUMENT HOLDER


- Keeps laparoscopic instruments close to operating field
- Reduces time required for locating instruments
- Multipocket sleeves to accommodate a variety of instrument sizes
- Available with tape to adhere to operating drape
- Available sterile and non-sterile

### STERILE


ITEM	DESCRIPTION	SIZE	PK
LK-700	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-800	4 Pockets + Tape	3.2 cm x 30.5 cm	25
LK-900	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25

### NON-STERILE


ITEM	DESCRIPTION	SIZE	PK
LK-700N	2 Pockets + Tape	6.5 cm x 30.5 cm	25
LK-900N	3 Pockets + Tape	3.2 cm x 30.5 cm (2), 6.5 cm x 30.5 cm	25



2 Pockets  
LK-700  
LK-700N



4 Pockets  
LK-800



3 Pockets  
LK-900  
LK-900N